

## Table of Contents

Chemical Reactions.....	2
Common Acids and Bases .....	7
Titration of Acids and Bases .....	12
Acidity in Soft Drinks .....	17
Properties of Solutions (Electrolytes) .....	23
Conductivity of Solutions .....	29
Action of a Buffer Solution.....	34
Freezing and Melting of Water .....	38
Intermolecular Forces and Evaporation .....	43
Water Cycle: Evaporation .....	47
Volume Changes and Gas Pressure .....	51
Gas Solubility .....	56
Combustion .....	60
Endothermic and Exothermic Reactions I.....	63
Endothermic and Exothermic Reactions II.....	67
Distillation Part I.....	71
Distillation Part II.....	75
Boyle’s Law .....	79
Producing Electricity .....	84
Solar Cells.....	88
Buffered Vitamin C Tablets .....	96
Study of Slow Release Food Supplements.....	101
Seed Respiration .....	105
Acid Rain .....	110

# Chemical Reactions

## Objective

- To learn about chemical reactions.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Oxygen Sensor: S98242-8ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- 250 ml beaker or cup
- Stirring rod
- Sodium sulfite (1/4 teaspoon)

## Discussion

A chemical reaction is a process that leads to the transformation of one set of chemical substances to another. Chemical reactions can be either spontaneous, requiring no input of energy, or non-spontaneous, often coming about only after the input of some type of energy, viz. heat, light or electricity. Chemical reactions encompass changes that strictly involve the motion of electrons in the forming and breaking of chemical bonds, although the general concept of a chemical reaction, in particular the notion of a chemical equation, is applicable to transformations of elementary particles, as well as nuclear reactions.

The substance/substances initially involved in a chemical reaction are called reactants. Chemical reactions are usually characterized by a chemical change, and they yield one or more products, which usually have properties different from the reactants.


Different chemical reactions are used in combination in chemical synthesis in order to get the desired product. In biochemistry, series of chemical reactions catalyzed by enzymes form metabolic pathways are performed within a cell.

Various NeuLog sensors can be used to monitor the changes in concentration of either reactants or products. For example, the Pressure Sensor can be used to monitor the

production and release of a gas in a chemical reaction. The pH Sensor can be used to follow an acid with a base (titration).

In this activity, you will use the Oxygen Sensor to follow a reaction between sodium sulfite<sup>1</sup> and oxygen dissolved in water.

The reaction is the following:



You will introduce the Oxygen Sensor into tap water, measure the oxygen concentration in % of dissolved oxygen and then add the sodium sulfite. You will then follow the changes in dissolved oxygen concentration.

## Procedure

### Experiment setup

1. Pour 200 ml tap water at room temperature into the beaker.
2. Place a paper towel under the sensor since some liquid may spill out. Remove the rubber protection from the Oxygen Sensor's  cap. Unscrew the cap, fill half of it with the included liquid (DO filling solution) and replace the cap.
3. Immerse the Oxygen Sensor in the beaker with water.



<sup>1</sup> The sodium sulfite is used in water treatment as an oxygen scavenger agent.

### Sensor setup

4. Connect the USB Module  to the USB port on your computer.
5. Connect the Oxygen Sensor  to the USB Module using a Data Cable.


### **Important:**

**You must wait at least 5 minutes after connecting the sensor to the USB Module before beginning calibration and measurements.**

6. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
7. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
8. This sets up the experiment parameters as follows:  
Experiment duration to 1 minute  
Sampling rate to 20 per second  
Sensor's mode: % in liquid
9. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
10. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 1 minute  
Set the sampling rate to 20 per second on the drop down menu.
11. Close the dialogue box.
12. Click on the 'Module Setup' icon  in the sensor's module box to open a dialogue box.



13. Since you are going to measure dissolved oxygen, set this mode by selecting the radio button next to % in liquid, as shown in the following picture:



14. Observe the measured values of dissolved oxygen in the module box. If 5 minutes have already passed since you connected the sensor to the USB module, the values will be stable. Whether you used the 'Load Configuration' function or not, click the 'Calibration' icon  to set the sensor to 100%<sup>2</sup>.
15. Close the module setup dialogue box.

### Testing and measurements

16. Click on the 'Run Experiment' icon  to start the measurement.
17. Allow the percentage of oxygen in water to stabilize (approximately 10 seconds).
18. Without stopping the measurement, add ¼ teaspoon sodium sulfite ( $\text{Na}_2\text{SO}_3$ ) to the water and stir continuously but carefully and far from the oxygen electrode.
19. Observe the changes in oxygen percentage.

## Summary Questions

- Describe what you see in the graph in general terms.
- Analyze and explain the changes in oxygen percentage.
- In terms of a mathematical function, which do you think describes better the behavior of the dissolved percentage of oxygen, linear, quadratic, exponential, etc?

<sup>2</sup> By doing this for simplicity, we are assuming 100% saturation of dissolved oxygen in this sample of water.

## Challenge

1. Use another amount of sodium sulfite.
2. Heat the water and repeat the experiment. How does this affect the rate of reaction?
3. Use the mathematical tools of NeuLog software to calculate the  $\ln$  of the measured concentration and see what kind of graph you receive.
4. What can you conclude about the chemical kinetics of the reaction, i.e. what is the order of this reaction?

# Common Acids and Bases

## Objective

- To determine the pH values of common household substances.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm – S98242-48ND
- USB Module - S98242-45ND 
- pH Sensor - S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Ring stand with utility clamp


## Materials

- Household solutions: tap water, vinegar, soft drink, ammonia, lemon juice, milk, soapy solution, baking soda solution, black coffee, etc.
- Paper cups
- Distilled or tap water

## Discussion

The pH measurement indicates if a liquid or solutions is acidic, basic or neutral. Many common household solutions contain acids and bases which play an essential role in our everyday lives. The pH scale is logarithmic and it shows the hydrogen ion  $[H^+]$  concentration in a solution, so acidity can be expressed using this scale that goes from 0 to 14. Neutral solutions have a pH 7, acidic solutions pH is less than 7 and basic solutions pH is more than 7.

In this experiment you will study the pH value of daily common solutions using a pH Sensor. You will determine the pH values of household substances such as lemon juice, vinegar, soda, black coffee, tap water, milk, soapy solution, baking soda.



## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Unscrew the cap of the storage solution and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
2. Put the pH Sensor in the tube holder on the stand and introduce it into the water stirring carefully to remove the soaking solution.



**Sensor setup**

3. Connect the USB Module  to the USB port on your computer.
4. Connect the pH Sensor  to the USB Module using a Data Cable.
5. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
6. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
7. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
8. This sets up a table for the results.
9. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
10. Click on the 'Experiment Setup' icon . This opens a dialogue box. Click the box next to table. Close the dialogue box.
11. Whether you used the 'Load Configuration' function or not, click on the Module Setup' icon  in the sensors module box to open a dialogue box.
12. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).



13. Close the Module Setup dialogue box.

**Testing and measurements**

14. Click on the 'Run Experiment' icon  to start the measurement.
15. Put the same amount of liquid, 20 ml, of each household solution into separate cups.
16. Introduce the pH Sensor into each solution, allow the pH to stabilize and then click the 'Single Step Mode' icon .
17. Remember to set aside a cup with tap water to rinse the sensor before inserting it into a new solution.
18. After each measurement, wait till the reading stabilized and write down the pH of each solution in the table.
19. At the end of the experiment, rinse the Ph sensor again, dip it in the soaking solution and screw the cap back in place.

**Data Collection:**

<b><i>Solution</i></b>	<b><i>pH</i></b>
Tap water	
Vinegar	
Soft drink	
Lemon juice	
Ammonia	
Milk	
Soapy solution	
Baking soda	
Black coffee	

## Summary Questions

1. Prepare a table of the solutions tested and record which are basic and which are acid.

## Challenge

1. Mix two or more acidic solutions and see how this affects the pH.
2. Mix one acidic solution with a basic one. See how this affects the pH of the mixture.

# Titration of Acids and Bases

## Objective

- To study the titration process and its applications.
- To follow changes in the pH during the titration process while adding a base to an acidic solution.
- To plot a graph of the titration process
  - The pH as a function of volume of the NaOH solution
- To determine the equivalence point of the titration by observing a graph.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- pH Sensor: S98242-24ND 

NOTE: Individual **NeuLog™** sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

These activity kits designed to explore titration using NeuLog sensors are available on our website:

Titration of Sodium Thiosulfate:	S02144
Titration of Polyprotic Acids:	S02145
Determination of Acid Concentration by Potentiometric Titration:	S02146

## Equipment and Accessories

- 50 ml beaker
- Ring stand with utility clamp
- A burette


## Materials

- 1 M of HCl
- 0.1 M of NaOH
- Tissue paper
- Water

## Discussion

Titration is an analytical tool to determine concentration of an acid or basic solution. It is based on the neutralization process. Neutralization occurs when the hydronium ion from an acid interacts with a hydroxide ion from a base, on a one-to-one basis, forming water in the process. A salt is always a byproduct of this type of reaction. Titration is the progressive addition of an acid to a base, or vice-versa, to achieve neutralization. The point at which the acid and base are in equivalent amounts is called the equivalence or end point.

One common example for acid-base titration is using a hydrochloric acid solution, HCl, with a basic sodium hydroxide solution, NaOH. In this experiment, you will titrate HCl using a NaOH solution with a known concentration and by using the end point; you will determine the concentration of HCl. In the reaction that takes place, hydrogen ions from the HCl react with hydroxide ions from the NaOH to produce water in a one-to-one ratio as shown in the following reaction:



At the beginning of the experiment the pH of the acidic solution is low. As we add the base the pH will gradually increase until it reaches the equivalence point. If we would add more of a base near the equivalence point, the pH would rapidly increase and then again would change gradually.

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble the following system:




2. Put the pH Sensor in the tube holder on the stand and introduce it in water stirring carefully to remove the soaking solution.
3. Pour 20 ml of water and 1 ml of 1M HCl into a clean beaker. Put the beaker in front of the stand.
4. Handle the acid and base with care, avoiding contact with the skin and clothes and wear personal protective products.
5. Prepare a diluted NaOH solution by mixing 1 ml 1M NaOH and 9 ml water in a separate disposable cup.
6. Attach the burette to the ring stand through the special burette holder. Fill the burette with the 0.1 M NaOH solution to the 0.0 ml mark.

### Sensor setup

7. Connect the USB Module  to the USB port on your computer.
8. Connect the pH sensor  to the USB Module using a Data Cable.
9. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
10. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
11. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
12. This sets up a table for the results.
13. Change the "Manual Values" column name to "Volume" in your table (clicking on the columns title will allow you to do so).
14. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.

15. Click on the 'Experiment Setup' icon . This opens a dialogue box: Click the box next to table. Change the "Manual Values" column name to "Volume" in your table (clicking on the columns title will allow you to do so).
16. Close the dialogue box.
17. Whether you used the 'Load Configuration' function or not, click on the 'Module Setup' icon  in the sensors module box to open a dialogue box.



18. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).
19. Close the module setup dialogue box.

### Testing and measurements

20. Take the pH Sensor out of the water and insert into the beaker.
21. Allow the pH to stabilize.
22. When the reading has stabilized, click the 'Single Step Mode' icon . Write the value 0 ml in the volume column in your table.
23. Begin the titration by adding 0.1 ml of the NaOH solution at a time using the burette.
24. Measure the pH each time clicking the 'Single Step Mode' icon .
25. Allow the pH to stabilize between measurements.

26. The new pH reading will appear in the table after each measurement; write the volume added next to each pH value.
27. When the pH value is about 3.5 begin adding one drop (instead of 0.1 ml). For each drop click the 'Single Step Mode'  icon to perform a new measurement. Write the volume next to the pH value that appears in the table.
28. When the pH is above 9, go back to adding 0.1 ml of the NaOH solution, and continue the pH measurements as before.
29. Keep adding the base until the pH value remains constant.

## Summary Questions

1. Plot a graph of the pH vs. volume of NaOH added in the titration. Click on the 'Experiment Setup' icon  and then click the 'XY Graph' tab. Select the 'Sensor' radio button. This will allow you to change the X in your graph to "Volume" using the down-arrow .
2. Follow the graph and find the volume of NaOH used in the titration. Also, find at what volume the largest jump in pH occurred. You can use the 'Show Cursors' icon . The values at the bottom left-hand corner of the graph indicate the X and Y for each cursor and the difference between X1 and X2 and between Y1 and Y2.
3. Calculate the volume of NaOH added to reach the equivalence point. To do so, find the volume before and after the jump, add them and divide by 2.

## Challenge

1. If you know the NaOH concentration and the volume of NaOH added to reach the end point, could you calculate the equivalent number of moles of NaOH?
2. Can you determine the concentration of the acid HCl in the beaker?
3. You can repeat the experiment with a different HCl concentration.

# Acidity in Soft Drinks

## Objective

- To monitor the change in pH, during titration, when adding a base to a soft drink acidic solution.
- To determine the equivalence point of the titration.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- pH Sensor: S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- 50 ml beaker
- Ring stand with utility clamp
- A Burette and a burette holder

## Materials

- 0.1 M NaOH
- Soft drinks: Coca Cola, Pepsi Cola, Sprite or 7-Up
- Tissue paper

## Discussion

The acidic content of many foods and beverages contributes significantly to their taste. Soft drinks often contain varying quantities of several acids, which give them their tart flavor. In cola products, these acids are predominantly carbonic (from carbonated water) and phosphoric acid. In sodas such as Squirt and 7-Up, the acids are carbonic and citric acid.

In this experiment, you will titrate different types of soft drinks with a basic sodium hydroxide solution, NaOH. The concentration of the NaOH solution is given and you will determine the

concentration of the different acids in the soft drinks. The point at which the acid and base are in equivalent amounts is called the end point. The following are examples of different types of acids undergoing a reaction with sodium hydroxide:




At the beginning of the experiment the pH of the acidic solution is low. As you add the base the pH will gradually increase until it reaches the equivalence point. If you add more base near the equivalence point the pH will rapidly increase and then again will change gradually.

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble the following system:




2. To prepare the soft drinks solution: Pour 20 ml of water into the beaker and add 20 ml of the soft drink you will test. Use each time a different sample. Allow the soft drink to stand until it loses the carbonic acid so that we do not have to take it into consideration.
3. In a separate cup, prepare a diluted NaOH solution by adding 5 ml of 1M NaOH to 50 ml of water.
4. Attach the burette to the ring stand with the special burette holder. Fill the burette with the 0.1 M NaOH solution to the 0.0 ml mark.

### Sensor setup

5. Connect the USB Module  to the USB port of your computer.
6. Connect the pH sensor  to the USB Module using a Data Cable.
7. Unscrew the cap of the storage solution and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
8. Put the pH Sensor in the tube holder on the stand and introduce it into the water stirring carefully to remove the soaking solution.
9. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
10. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
11. This sets up a table for the results.
12. Change the "Manual Values" column name to "Volume" in your table (clicking on the columns title will allow you to do so).
13. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
14. Click on the 'Experiment Setup' icon . This opens a dialogue box: Click the box next to table.
15. Change the "Manual Values" column name to "Volume" in your table (clicking on the columns title will allow you to do so).

16. Close the dialogue box.
17. Whether you used the 'Load Configuration' function or not, click on the 'Module Setup' icon  in the sensors module box to open a dialogue box.



18. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).
19. Close the module setup dialogue box.

### Testing and measurements

20. Take the sample you prepared in step 2 of the experiment setup.
21. Insert the pH Sensor into the beaker. Allow the pH to stabilize.
22. When the reading has stabilized, click the 'Single Step Mode'  icon. Write the value 0 ml in the volume column in your table.
23. Begin the titration by adding 0.5 ml of the NaOH solution at a time using the burette.
24. Measure the pH each time clicking the 'Single Step Mode'  icon.
25. Allow the pH to stabilize between measurements.
26. The new pH reading will appear in the table after each measurement; write the volume added next to each pH value.
27. When the pH value is about 5, begin adding 0.1 ml (instead of 0.5 ml). For each drop, click the 'Single Step Mode'  icon to perform a new measurement. Write the volume next to the pH value that appears in the table.

28. Continue until the pH is above 9. Save the data at the end of the experiment.
29. Rinse the pH Sensor with water, without disconnecting it from the USB Module.
30. Rinse the beaker and prepare another sample.
31. Repeat steps 19-28 for the new sample.
32. Repeat the experiment with all the samples prepared.
33. At the end of the experiments, rinse the pH Sensor and return it to the storage solution.

**Data Table:**

Sample	volume of NaOH added (ml)	Moles of NaOH added	Moles of the Acid

## Summary Questions

For each one of the soft drinks you tested, calculate the volume of NaOH added to reach the equivalence point and write it in your data table:

1. Plot a graph of the pH vs. volume of NaOH added in the titration. Click on the 'Experiment Setup' icon  and then click the **XY Graph** tab. Select the **Sensor** radio button. This will allow you to change the X in your graph to "Volume" using the down-arrow .
2. Follow the graph and find the volume of NaOH used in the titration. Also, find at what volume the largest jump in pH occurred. You can use the 'Show Cursors' icon . The values at the bottom left-hand corner of the graph indicate the X and Y for each cursor and the difference between X1 and X2 and between Y1 and Y2.
3. Calculate the volume of NaOH added to reach the equivalence point. To do so find the volume before and after the jump, add them and divide by 2.

4. Calculate the number of moles of NaOH added to reach the end point for each one of the soft drinks you tested. Write it in your data table.

## Challenge

1. Could you determine the moles of acid in each drink?
2. Could you determine the concentration of the acid in each drink?

## Properties of Solutions (Electrolytes)

### Objective

- To measure the conductivity of common solutions.
- To write equations for the dissociation of compounds in water.
- To study the relationship between molecular dissociation products and their conductivity.

### Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Conductivity Sensor: S98242-5ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

### Equipment and Accessories

- Ring stand
- Utility clamp

### Materials

- Disposable cups
- Disposable tea spoons
- Distilled water
- 0.05 M of the following compounds: NaCl, CaCl<sub>2</sub>, HCl, C<sub>2</sub>H<sub>5</sub>OH (ethanol) and sugar
- Tissue paper

### Discussion

You've probably been warned about not standing in water during lightning, or keeping electrical appliances away from the sink or tub, but did you know that the electrical conductivity of *pure* water is not good? The materials dissolved in water are those which conduct electricity. Most bodies of water - even tap water - contain electrolytes, compounds that ionize in water, and this produces a solution that conducts electricity.

The ability of a solution to conduct electricity (its *conductivity*) is measured by a device called conductivity meter. Conductivity readings can be used to estimate how many electrolytes are dissolved in the water, and hence tell how clean a sample of "purified" water is. Conductivity readings are used by water treatment plants to grade water quality (called the "total dissolved solids") to determine the water's suitability for human consumption.

Electrolytic substances are classified as strong or weak, according to how readily they dissociate into conducting ions. The stronger the electrolyte is in the substance, the higher the conductivity readings.

In this experiment, you will discover some properties of strong electrolytes, weak electrolytes and non-electrolytes by observing the behavior of substances in aqueous solutions. You will observe several factors that determine whether a solution will conduct electricity and the relative magnitude of conductivity.

The characteristics of the three types of electrolytes are:

Strong electrolyte – Usually ionic compounds: a strong electrolyte dissociates or ionizes completely or almost completely to form free mobile ions in the solution or molten form. The greater the availability of free mobile ions in an electrolyte, the greater its capacity to carry or conduct a current.

Example:  $\text{NaCl} \rightarrow \text{Na}^+ + \text{Cl}^-$

Weak electrolyte – dissolves partially, most molecular acids, (some molecular acids are strong).

Non electrolyte – Most molecular compounds: They do not dissociate to form ions and therefore do not conduct electricity.


For example, ethanol:  $\text{C}_2\text{H}_5\text{OH}(\text{l}) \rightarrow \text{C}_2\text{H}_5\text{OH}(\text{aq})$ .

In this experiment you will study different solutions to determine if they conduct electricity or not and what type of electrolyte they are.

## Procedure

### Experiment setup

1. Hold the sensor using the ring stand and clamp as in the picture:



2. Prepare the solutions as follows:
  - a. Dissolve 0.5 gram (a quarter of a tea spoon) NaCl salt in 10 ml distilled water; this solution is 1 M. Take 1 ml of the 1M solution and add 20 ml distilled water to produce the 0.05 M solution needed.
  - b. Dissolve 1 gram (half tea spoon) of CaCl<sub>2</sub> salt in 10 ml distilled water; this solution is 1 M. Take 1 ml of the 1M solution and add 20 ml distilled water to produce the 0.05 M solution needed.
  - c. Dissolve 1 ml 1 M HCl in 20 ml distilled water to produce the 0.05 M solution needed.
  - d. Dissolve 1 ml Ethanol (concentrated) in 20 ml distilled water to produce the 0.05 M solution needed.
  - e. Dissolve 3.5 grams (two tea spoons of sugar) in 10 ml distilled water; this solution is 1 M.
  - f. Take 1 ml of the 1M solution and add 20 ml distilled water to produce the 0.05 M solution needed.

### Sensor setup

3. Connect the USB Module  to the USB port on your computer.
4. Connect the 'Conductivity' Sensor  to the USB Module using a Data Cable.
5. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
6. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
7. This opens a data table; change the "Manual Values" column name to "Sample" (clicking on the title Manual Values will allow you to change it). It also sets the sensor's units to  $\mu\text{s}/\text{cm}$ .
8. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
9. Click on the 'Experiment Setup' icon . This opens a dialogue box.
10. Select the table display. Change the "Manual Values" column name to "Sample" (clicking on the title Manual Values will allow you to change it).
11. Close the dialogue box.
12. Click on the 'Module Setup' icon  to open a dialogue box. Set the sensor's units by selecting the radio button next to  $\mu\text{s}/\text{cm}$ .
13. Close the Module Setup dialogue box.

### Testing and Measurements

14. Introduce the 'Conductivity' sensor  in the cup with distilled water. Make sure the sensor is completely submerged in the solution to be tested. The reading in the sensor's module box on the left side of the screen should indicate no conductivity (low reading, close to zero).
15. After each measurement, wash the Conductivity Sensor carefully with water.

16. Measure the conductivity of distilled water. When the reading has stabilized, click the

‘Single Step Mode’  icon. Write the sample name in your table.

17. Measure the conductivity for each solution as follows:

- Stir the cup contents.
- Dip the sensor and wait for the conductivity reading to stabilize;
- Click the ‘Single Step Mode’  icon to measure and record the sample name in your data table.

18. Between solutions, remember to wash the sensor with distilled water and dry it using tissue paper.

19. Save the data at the end of the experiment.

**Data Collection:**

<i><b>Solution</b></i>	<i><b>Conductivity (<math>\mu\text{S}/\text{cm}</math>)</b></i>
<b>Distilled water</b>	
<b>Tap water</b>	
<b>HCl</b>	
<b>C<sub>2</sub>H<sub>5</sub>OH</b>	
<b>Sugar</b>	
<b>NaCl</b>	
<b>CaCl<sub>2</sub></b>	

## Summary Questions

1. Observe the conductivity values measured and divide the electrolytes in groups.
2. Write an equation for the dissociation of each compound.
3. Based on the dissociation, could you explain the difference in conductivity?

## Challenge

1. How do you explain the difference in the conductivity of tap water from different countries or even different cities in the same country?
2. Test water from different sources (tap water, filtered water, mineral water, distilled water) and compare the conductivity measurements.

# Conductivity of Solutions

## Objective

- To study the relationship between conductivity and concentration of solutions.
- To study the relationship between the number of ions produced in dissociation and the conductivity.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Conductivity Sensor: S98242-5ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Ring stand
- Utility clamp

## Materials


- 100 ml cup
- Distilled water
- 1 M of the following compounds: NaCl and CaCl<sub>2</sub>
- Tissue

## Discussion

The ability of a solution to conduct electricity (its *conductivity*) is measured by a device called a conductivity meter. The conductivity readings can be used to estimate how many electrolytes are dissolved in water, and hence tell how clean a sample of "purified" water is. Conductivity readings are used by water treatment plants to grade water quality (called the "total dissolved solids") to determine the water's suitability for human consumption.

Electrolytic substances are classified as strong or weak according to how readily they dissociate into conducting ions. The stronger the electrolyte is, the higher the conductivity.

When solid sodium chloride dissolves in water it releases ions as we see in the following equation:




In this experiment you will gradually increase the concentration of an ionic compound and measure its conductivity. First, you will study the effect of an increased concentration on conductivity by gradually adding concentrated NaCl. Then, you will study the effect of adding the same concentration of another solution, with a different number of ions as they dissolve in water, such as: Calcium Chloride,  $\text{CaCl}_2$ .

## Procedure

### Experiment setup

1. Hold the sensor using the ring stand and a clamp as in the picture.



2. Add 70 ml of distilled water to a clean 100-ml cup and put it in front of the sensor.

3. Prepare the solutions as follows:
  - a. Dissolve half a tea spoon of NaCl salt, about 1 gram, in 20 ml of distilled water for 1 M NaCl
  - b. Dissolve one tea spoon of CaCl<sub>2</sub> salt, about 2 grams, in 20 ml distilled water for 1 M CaCl<sub>2</sub>

### Sensor setup

4. Connect the USB Module  to the USB port on your computer.
5. Connect the Conductivity Sensor  to the USB Module using a Data Cable.
6. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
7. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
8. This opens a data table; change the "Manual Values" column name to "Sample" (clicking on the title Manual Values will allow you to change it). It also sets the sensor's units to **µs/cm**.
9. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
10. Click on the 'Experiment Setup' icon . This opens a dialogue box.
11. Select the table display. Change the "Manual Values" column name to "Sample" (clicking on the title Manual Values will allow you to change it).
12. Close the dialogue box.
13. Click on the 'Module Setup' icon  to open a dialogue box. Set the sensor's units by selecting the radio button next to **µs/cm**.
14. Close the Module Setup dialogue box.

### Testing and measurements

15. Introduce the Conductivity Sensor  in the cup with distilled water. Make sure the sensor is completely submerged in the solution to be tested. The reading in the

sensor's module box on the left side of the screen should indicate no conductivity (low reading, close to zero).

16. Prepare a cup with clean water for washing the Conductivity Sensor carefully after each measurement.
17. Measure the conductivity of distilled water. When the reading has stabilized, click the 'Single Step Mode'  icon. Write the value of conductivity as reference 0 ions added.
18. Add the salt solution one drop at a time, measuring the conductivity for each added drop as follows:
  - a. Add one drop of 1M NaCl to the distilled water. Stir the cup contents. Deep the sensor and once the conductivity reading has stabilized, click the 'Single Step Mode'  icon.
  - b. Record the value in your data table.
  - c. Continue drop by drop until a total of 7 drops has been added.
19. Wash the sensor with water end dry it with tissue paper.
20. Save the data.
21. Plot a graph of concentration vs. volume. Click on the 'Experiment Setup' icon  then click on the **XY Graph** tab if you chose the "Sensor" option, you will be able to change the X in your graph to "Volume" by pressing the down-arrow  .
22. Calculate the slop value. Press the 'Show Functions'  icon;
  - a. Click on the functions tab and select Linear Fit from the second drop down menu.
  - b. Click the 'Calculate Function'  icon to perform the calculation.
  - c. Observe, in the box, the formula with the best fit line.
  - d. Record the value of the slope in your data table.
  - e. Erase the data before the next test by clicking the 'Clear Experiment Results'  icon.
23. Repeat Steps 14-22, using this time 1.0 M CaCl<sub>2</sub> solution in step 18.

**Data Collection:**

<i>Solution Volume in # drops</i>	<i>Conductivity (μS/cm) NaCl</i>	<i>Conductivity (μS/cm) CaCl<sub>2</sub></i>
0		
1		
2		
3		
4		
5		
6		
7		

<i>Solution</i>	<i>Slope</i>
1 M NaCl	
1 M CaCl <sub>2</sub>	

### Summary Questions

1. What are the dissociation products of NaCl and CaCl<sub>2</sub> in water? Write the chemical dissociation equation for both salts.
2. Observe the conductivity graph as a function of concentration of the NaCl solution. What kind of mathematical relationship exists between conductivity and concentration?
3. Does the same relationship exist between conductivity and concentration for the CaCl<sub>2</sub> solution? Which solution has the higher slope and why?

### Challenge

1. If the conductivity of NaCl is known, could you predict its concentration?

# Action of a Buffer Solution

## Objective

- To investigate the action of a buffer solution.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- pH Sensor: S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- 20 ml beaker
- Droppers

## Materials

- Alka Seltzer – buffer
- 1 M NaOH and HCl solutions  
(household materials such as vinegar and ammonia could be used instead)

## Discussion

In chemistry, a buffer is a system, usually an aqueous (water) solution, which resists having its pH changed when an acid or a base is added to it.

Buffer solutions are used as means of keeping pH at a nearly constant value in a wide variety of chemical applications. Many life forms thrive only in a relatively small pH range. An example of a buffer solution is blood.

A Buffer Solution contains a weak acid and a salt of the same acid or a weak base and a salt of the same base and has the property of reacting both with acids and bases.


Through this activity, we shall discover the main characteristics of such mixtures. The first part of the experiment will allow us to study how water reacts to the addition of an acid or a base. In the second part, the experiment will be repeated using this time an aqueous solution with a buffer (Alka Seltzer).

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble a system as in the picture below.




2. Pour 15 ml of water into a beaker.
3. Take out the pH Sensor and unscrew the cap of the storage solution. Take out the Sensor and set the solution aside. Raise the cap to the top of the Sensor.
4. Insert the Sensor into the water and stir carefully to remove the storage solution.

### Sensor setup

5. Connect the USB Module  to the USB port on your computer.

6. Connect the pH Sensor  to the USB Module using the Data Cable.
7. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
8. Observe the sensor's Module Box on the left side of the screen. Allow the pH value to stabilize.
9. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
10. This sets up the experiment parameters as follows:  
Experiment duration to 3 minutes  
Sampling rate to 10 per second
11. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
12. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 3 minutes  
Set the sampling rate to 10 per second on the drop down menu.
13. Close the dialogue box.
14. Whether you used the Load Configuration function or not, click on the 'Module Setup' icon  in the sensors Module Box to open a dialogue box.
15. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).
16. Close the module setup dialogue box.



## Testing and measurements

Note: In order to understand the effect of Alka-zeltzer, we will first measure the pH variation while adding an acid to water, then a base. Next, we will repeat the procedure adding the acid or the base to water with an Alka-zeltzer tablet dissolved.

17. Click on the 'Run Experiment' icon  to start the measurement.
18. Without stopping the measurement, use a dropper to add about 1 ml hydrochloric acid (HCl) to the water in the beaker and stir slowly. Wait until the end of the measurement.
19. Click on the 'Freeze Current Graphs' icon  .
20. Wash the beaker and electrode with water. Pour 15 ml of fresh water into the clean beaker and insert the pH Sensor.
21. Click on the 'Run Experiment' icon  to start a new measurement.
22. Allow the pH to stabilize.
23. Without stopping the measurement, use another dropper to add about 1 ml sodium hydroxide (NaOH) to the water and stir slowly. Wait until the end of the measurement.
24. Freeze the graph as in step 17.
25. Observe the changes in the pH in both cases.
26. Carry out the same two experiments, using water with a buffer, i.e. In each experiment, add half an Alka Seltzer tablet to the water. Wait until it dissolves well and repeat steps 16-22.

## Summary Questions

1. Observe the four graphs you obtained and describe in detail what happened to the pH in each case.
2. Investigate what are the two chemical components needed to produce a buffer.
3. Describe in chemical terms how a buffer works.

## Challenge

1. Investigate what is the difference between a buffer medicine and an antacid.
2. Investigate what is the importance of buffer solutions.

# Freezing and Melting of Water

## Objective

- To understand phase changes.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Ring stand
- Utility clamp
- Test tube 20 x 150 mm
- 400 mL plastic beaker
- 10 mL graduated cylinder
- Stirring rod

## Materials

- 300 mL crushed ice
- 200 g salt

## Discussion



In physical science, freezing or solidification is the process in which a liquid turns into a solid when cold enough. The freezing point is the temperature at which this happens. Melting, the process of turning a solid into a liquid is almost the exact opposite of freezing.

All known liquids undergo freezing when the temperature is lowered enough, with the sole exception of liquid helium, which remains liquid at absolute zero and can only be solidified under pressure. For most substances, the melting and freezing points are the same temperature; however, certain substances possess differing solid-liquid transition

temperatures. For example, agar melts at 85 °C (185 °F) and solidifies from 31 °C to 40 °C (89.6 °F to 104 °F); this process is known as hysteresis.

In this experiment you will use a sensor to follow the temperature changes of liquid water starting at room temperature and cooling it up to temperatures of -10 °C. Then you will do a second measurement starting at -10 °C and letting the sample heat up to room temperature again. This will allow you to get both the cooling and the heating curve of water passing through its solid to liquid phase changes.

Use the following graphs to draw a hypothesis of how you think the temperature measurement will look as a function of time. Do it both for the first part of the experiment (freezing of water) and for the second part (melting of water).



## Procedure

### Sensor setup

1. Connect the USB Module  to the USB port on your computer.
2. Connect the Temperature Sensor  to the USB Module using a Data Cable.
3. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
4. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .

5. This sets up the experiment parameters as follows:  
 Experiment duration to 15 minutes  
 Sampling rate to 10 per second  
 Ymax =25, Ymin = -25, Y axis position = -7
6. This icon will setup the experiment parameters so you can skip the '**Module and experiment setup**' steps, as follows:  
 Experiment duration to 15 minutes.  
 Sampling rate to 10 per second
7. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
8. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
 Set the Experiment Duration to 10 minutes  
 Set the sampling rate to 60 per minute on the drop down menu.
9. Close the dialogue box.
10. Click the 'Sensor Setup'  icon in the sensors module box on the left side of the screen to open a dialogue box.
11. Select the graph tab in the opened dialogue box.
12. Use the arrows to set the Y max. = 25, then the Y min = -25 and the Y axis position = -7.



13. Close the module setup dialogue box.

### Experiment setup

14. Assemble a system as in the picture below:



15. Fill the beaker with 4 parts of crushed ice and 1 part of salt (up to 3/4 of the beaker). Stir well. Continue to stir from time to time during the experiment.
16. Pour 5 ml of water into the test tube; fasten the test tube to a ring stand with the utility clamp. **The test tube should be above the plastic beaker.**
17. Insert the temperature sensor into the test tube with water.

### Testing and measurements

#### Freezing of water

18. Click on the 'Run Experiment'  icon to start the measurement.
19. Lower the test tube into the crushed ice in the beaker. Verify that all the water in the test tube is immersed in the ice.
20. Keep stirring from time to time the ice / salt mixture.
21. Also, carefully stir the water with the sensor a couple of times during the first minutes.
22. The measurement will end after 15 minutes. Observe the graph.

23. Save your data.

#### Melting of water



24. Click on the 'Run Experiment' icon to start a new measurement.

25. Raise the test tube and fasten it in a position above the plastic beaker.

26. Throw away the ice / salt mixture and fill the beaker with warm water.

27. **After about 10 minutes**, lower the test tube into the warm water bath.

28. When the measurement is over, observe the resultant graph.

29. Save the data.

### Summary Questions

1. Analyze both graphs and determine which are the freezing and melting points of water.
2. You may see in your resulting plots that the slope of cooling (or heating) is not the same. Divide the graph in three characteristically different zones and explain what each zone in the graph represents.

### Challenge

1. Design and perform a similar experiment with other chemicals, such as tert-Butanol or Lauric acid. (Investigate their melting points).

# Intermolecular Forces and Evaporation

## Objective

- To study the changes in temperature caused by the evaporation of several solvents.
- To relate the temperature difference due to evaporation to the strength of intermolecular attraction.

## Modules and Sensors

- Computer with NeuLog™ Software
- 3 Data Cables, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- 3 Temperature Sensors: S98242-31ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Materials

- Two small rubber bands
- Cotton pads
- Alcohol
- Acetone

## Discussion

During evaporation, the magnitude of a temperature decrease is an indicator of the force of attraction between molecules of a liquid. The greater the temperature drop, the faster the liquid evaporates and the weaker its intermolecular attraction is.

When changing from liquid to gaseous state, there is an increase in the potential energy of molecules. As a result, there is a decrease in the average kinetic energy of the molecules left behind in the liquid and thus, there is a temperature decrease. This phenomenon can be felt when stepping out of the shower in the morning since there is a slight drop in the temperature of the skin. The evaporation of water from the skin absorbs energy from the body. In this activity we will measure the change in temperature for three types of compounds and link them to perspiration.

## Procedure

### Experiment setup




1. Wrap the tips of two temperature sensors with a small piece of cotton and secure it with a small rubber band. Wrap the cotton evenly.

### Sensors setup

2. Connect the USB Module  to the USB port on your computer.
3. Connect one Temperature Sensor  to the USB Module using a Data Cable.
4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
5. Click on the 'Tools' icon  in the main icons bar to reveal the Tools sub icon bar.

- Bring the cursor to 'Set sensor ID Number' icon . Set the number to 2 and click the icon. A new sensors search will start automatically.
- Your sensor should appear now with ID = 2 (you will see the number 2 in the upper left part of the sensor's module box).



- Disconnect the sensor from the USB Module and connect another Temperature Sensor. Repeat step 4 with the new sensor assigning the ID number as 3 this time.
- Leave this sensor connected to the USB Module and connect the other two sensors forming a chain (any free socket will do).
- Click the 'Search for Sensors' icon  again. You should see the three connected sensors, each with a different ID number: 1, 2 and 3.
- If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click on the 'Load Configuration' icon .
- This sets up the experiment parameters as follows:  
Experiment duration to 2 minutes  
Sampling rate to 10 per second
- If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.

- Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 2 minutes  
Set the sampling rate to 10 per second on the drop down menu.

- Close the dialogue box.

### **Testing and measurements**

- Insert the tip of sensor with ID 1 in the alcohol carefully so that the cup does not tip over.
- Insert the tip of sensor with ID 2 in the acetone.

18. The third sensor is for control.
19. Leave the sensors in the liquids for about ten seconds.
20. Click on the 'Run Experiment' icon  to start the measurement.
21. Wait about 10 seconds and then remove the sensors simultaneously from the liquids.
22. Put the three sensors on the table with their tips coming out (see picture in previous page).
23. Observe the graph.

## Summary Questions

1. Calculate the temperature difference for each of the liquids.

## Challenge

1. Could you explain the different changes of evaporation temperature for the liquids tested in terms of intermolecular forces of attraction?

# The Water Cycle: Evaporation

## Objective

- To study changes in temperature caused by the evaporation of solvents.
- To relate the temperature difference due to evaporation.

## Modules and Sensors

- Computer with NeuLog™ Software
- 3 Data Cables, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- 3 Temperature Sensors: S98242-31ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Materials

- Two small rubber bands
- Cotton pads
- Alcohol
- Acetone

## Discussion

Evaporation is the process by which water changes from a liquid to a gas or vapor. Evaporation is the primary pathway that water moves from the liquid state back into the water cycle as atmospheric water vapor. Studies have shown that the oceans, seas, lakes, and rivers provide nearly 90 percent of the moisture in the atmosphere via evaporation, with the remaining 10 percent being contributed by plant transpiration.

Heat (energy) is necessary for evaporation to occur. Energy is used to break the bonds that hold water molecules together, which is why water easily evaporates at the boiling point (212° F, 100° C) but evaporates much more slowly at the freezing point. Net evaporation occurs when the rate of evaporation exceeds the rate of condensation. A state of saturation exists when these two process rates are equal, at which point the relative humidity of the air is 100 percent.

The process of evaporation removes heat from the environment, which is why water evaporating from your skin cools you. Also, your arm gets very cold when a physician rubs it with alcohol. In this experiment you will use a temperature sensor to study the changes in temperature during evaporation of alcohol or acetone instead of water since the evaporation process occurs faster.

## Procedure


### Experiment setup



1. Wrap the tips of two temperature sensors with a small piece of cotton and secure it with a small rubber band. Wrap the cotton evenly.

### Sensor setup

2. Connect the USB Module  to the USB port on your computer.
3. Connect one Temperature Sensor  to the USB Module using a Data Cable.
4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .

5. Click on the 'Tools' icon  in the main icons bar to reveal the Tools sub icon bar.
  - a. Bring the cursor to 'Set Sensor ID Number' icon .
  - b. Set the number to 2 and click the icon.
  - c. A new sensors search will start automatically;
  - d. your sensor should appear now with ID = 2
  - e. you will see the number 2 in the upper left part of the sensor's module box
6. Disconnect the sensor from the USB Module and connect another Temperature Sensor. Repeat step 4 with the new sensor assigning the ID number as 3 this time.
7. Leave this sensor connected to the USB Module and connect the other two sensors forming a chain (any free socket will do).
8. Click the 'Search for Sensors' icon  again. You should see the three connected sensors, each with a different ID number: 1, 2 and 3.
9. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
10. This sets up the experiment parameters as follows:  
 Experiment duration to 2 minutes  
 Sampling rate to 10 per second
11. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
12. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
 Set the Experiment Duration to 2 minutes  
 Set the sampling rate to 10 per second on the drop down menu.
13. Close the dialogue box.

### **Testing and measurements**

14. Insert the tip of sensor with ID 1 in the alcohol carefully so that the cup does not tip over.
15. Insert the tip of sensor with ID 2 in the acetone.
16. The third sensor is for control.

17. Leave the sensors in the liquids for about ten seconds.
18. Click on the 'Run Experiment' icon  to start the measurement.
19. Wait about 10 seconds and then remove the sensors simultaneously from the liquids.
20. Put the three sensors on the table with their tips coming out (see picture in previous page).
21. Observe the graph.

## Summary Questions

1. What happened to the temperature, did it increase, decrease or stayed stable?
2. Calculate the temperature difference.
3. Compare your graphs.
4. How could we relate the results to evaporation?

## Challenge

1. Could you explain the different changes of evaporation temperature for the liquids tested in terms of intermolecular forces of attraction?

# Volume Changes and Gas Pressure

## Objective

- To study the behavior of gases.
- Use a Pressure Sensor to measure the change in pressure as a function of change in volume.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Pressure Sensor: S98242-13ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Gas Law Device Kit: S02107  
(Volume of gas in the device is 0.055 L)

## Discussion

An ideal gas can be characterized by three parameters: volume, pressure and temperature, for a certain amount of gas. The law that presents the relationship between volume, pressure and temperature is called the Ideal Gas Law:  $PV = nRT$  where  $nR$  is constant when mass of the gas is constant; the temperature is in Kelvin degrees.

The law says that volume is inversely proportionate to pressure:  $PV = k$  where  $k$  is constant if the temperature and mass of the gas are constant.

When pressure increases, the volume decreases; if pressure decreases the volume increases. It is not necessary to know the exact value of the constant  $k$  in order to use the law; if the amount of gas and its temperature are kept constant, the relationship:  $P_1V_1 = P_2V_2$  must be fulfilled.

In this activity, we will change the volume of a gas container in order to measure how it affects its pressure.

## Procedure

### Experiment setup

1. Assemble a system as in the picture below.



2. Take out the kit for the study of gases.
3. Pull the syringe's plunger up to the 60 ml position and insert the syringe into one of the kit's perforated caps.

### Sensor setup

4. Connect the USB Module  to the USB port on your computer.
5. Connect the Pressure Sensor  to the USB Module using the Data Cable.
6. Insert the Pressure Sensor's inlet connector into one of the kit's perforated caps, as in the picture.

7. Make sure that the third cap is sealed and verify that there is no air leak by pressing the syringe's plunger and releasing it so that it goes back up. If necessary, take out the syringe, bring the plunger to the 60 ml mark and put it back in place.
8. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
9. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
10. This sets up the experiment parameters as follows: Sensor's units: Atm
11. If you did not use the 'Load Configuration' function , click on the 'Module Setup' icon  in the sensor's module box to open a dialogue box. Set the sensor's units by selecting the radio button next to Atm.

Your screen should look like the picture below:



### Testing and measurements

12. Click the 'On Line Experiment' icon  in the main icon bar. Click on the 'Experiment Setup' icon . This opens a dialogue box. Observe that the graph display is selected by default. Select the table display.
13. Click the 'Single Step' icon  to measure the gas pressure value in the container + the syringe (atmospheric pressure in your surroundings). Observe the graph and the value on the table.




14. Insert the syringe's plunger to the 50 ml position and click the 'Single Step' icon again for a new measurement.
15. Continue to decrease the volume to 40, 30 and 20 ml.
16. Add the volume values to the table as follows:
  - a. Click on Manual values to change this title;
  - b. Type Volume [L] (It is very important to write the units inside square parenthesis since the program recognizes and uses them for math calculations which affect the units.).
  - c. Press enter.
  - d. Type 0.115 in the first space of the volume column of your table (55 ml gas volume in the device + 60 ml gas volume in the syringe = 115 ml = 0.115 L).
  - e. Continue writing the values: 0.105, 0.095, 0.085, 0.075.

The next steps indicate how to calculate the reciprocal value of the volume.

17. Click the 'Functions' icon  and then click the math tabulator in the opened dialogue box.
18. Go to the A drop down menu and choose the Volume option.
19. Choose Function 1/A from the third drop down menu.
20. Click 'Calculate Math' icon .
21. Observe the table and see that a new column was created, titled 1/(Volume) [1/L] and under it the calculated values.
22. Observe also the graph with calculated values as a function of samples.

Your screen should look like the one below:



The next steps indicate how to plot a graph of pressure as a function of 1/(Volume).

23. Click the 'Experiment Setup' icon  and access the tabulator Graph XY.
24. Select the sensor radio button.
25. Select 1/Volume option from the drop down menu. You now have a graph of Pressure as a function of 1/Volume.

Our last step will be to calculate the linear fit for the points we received in the graph.

26. Click again the 'Functions' icon . Verify that the dialogue box opens in the functions tabulator and also that the functions in the second drop down menu is linear fit.
27. Click the 'Calculate Function' icon  in order to calculate the linear fit. Observe the graph with the correspondent line and in the functions window, the linear function calculated.
28. You can use the option of extrapolation to zero by clicking this  icon. Does the line cross the zero value? If not, why?

## Summary Questions

1. Use other volumes and compare your results.

# Gas Solubility

## Objective

- To learn about gas solubility.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Oxygen Sensor: S98242-8ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- A ring stand with utility clamp
- A 50 ml beaker or a cup
- A syringe

## Discussion

Solubility is a physical property that describes the ability of a given substance, the solute, to dissolve into another substance, the solvent. Solubility is measured in terms of the maximum amount of solute dissolved in a solvent at equilibrium. The resulting solution is called saturated.

The solubility of a gas in a liquid depends on temperature, the partial pressure of the gas over the liquid, the nature of the solvent and the nature of the gas. The most common solvent is water.

The concentration of dissolved gas depends on the partial pressure of that gas. The partial pressure controls the number of gas molecule collisions with the surface of the solution. If the partial pressure is doubled, the number of collisions with the surface will double. The increased number of collisions produces more dissolved gas.

The solubility of a gas such as oxygen, in a liquid, decreases when the temperature rises. The existing relationship between temperature and solubility is proportionally inverse; that is, the lower the temperature of the liquid, the more the gas will dissolve. However, other gases in

other solvents do not necessarily become less soluble at high temperatures.  $H_2$ ,  $N_2$ ,  $CO$ ,  $He$ , and  $Ne$  solubility actually rises with a higher temperature, in common organic solvents. In this activity we will saturate water with Oxygen (g).

## Procedure

### Experiment setup

1. Pour 20 ml tap water at room temperature into the beaker.
2. Remove the rubber protection from the Oxygen Sensor's cap. Unscrew the cap, fill half of it with the included liquid and close the cap.
3. Attach the Oxygen Sensor  to the ring stand.




### Sensor setup

4. Connect the USB Module  to the USB port on your computer.
5. Connect the Oxygen Sensor  to the USB Module using the Data Cable.

### Note:

**You must wait at least 5 minutes after connecting the sensor to the USB Module before beginning calibration and measurements.**

6. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
7. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
8. This sets up the experiment parameters as follows:  
 Experiment duration to 7 minutes  
 Sampling rate to 60 per minute  
 Sensor mode: % in air.
9. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
10. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
 Set the Experiment Duration to 7 minutes  
 Set the sampling rate to 60 per minute on the drop down menu.
11. Close the dialogue box.
12. Whether you used the 'Load Configuration' function or not, you must prepare the dissolved oxygen sensor for measurement:
  - a. Carefully remove the rubber protection from the sensor's cap.
  - b. Unscrew the cap, fill half of it with the included liquid (DO filling solution) and close it.
  - c. The calibration of the Oxygen Sensor in % in liquid assumes the saturation is 100%.
  - d. In order to reach this condition, introduce the tip of the Oxygen Sensor in a container with tap water and inject air with the 50 mL syringe.
  - e. Pump about five times.
13. Immediately afterwards, look for the Oxygen Sensor's module box on the left side of the screen and click on the 'Module Setup' icon  to open a dialogue box.



14. Verify the mode is "% in liquid". Click the 'Calibration' icon  which sets the oxygen value to 100% (of oxygen saturation in water).
15. After pumping air, the reading of dissolved oxygen decreases which is normal. The reading shows the real value of dissolved oxygen in the experimental conditions. The oxygen leaves the water and goes to the air.
16. Close the module setup dialogue box.

### Testing and measurements

17. Click on the 'Run Experiment' icon  to start the measurement.
18. Allow the percentage of oxygen in water to stabilize (approximately 30 seconds).
19. Without stopping the measurement, insert the syringe and pump air into your sample.
20. Keep on pumping air for about 2-3 minutes.
21. Without stopping the measurement, take out the syringe from your sample.
22. Observe the changes in the oxygen percentage.

### Summary Questions

1. Describe what you see in the graph.
2. Analyze and explain the changes in the oxygen percentage
3. Why when you stop pumping the concentration of dissolved oxygen decreases?

### Challenge

1. Explain, giving three examples, the importance of aqueous solutions.
2. What is the influence of temperature and pressure in the solubility of a gas in water?

# Combustion

## Objective

- To determine the factors that influences the chemical process of combustion.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Oxygen Sensor: S98242-8ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Gas Law Device Kit: S02107  
(Volume of gas in the device is 0.055 L)


## Materials

- A candle

## Discussion

**Combustion** or **burning** is a complex sequence of exothermic chemical reactions between a fuel and an oxidant accompanied by the production of heat or both heat and light in the form of either glow or flames, appearance of light flickering.

Many organic compounds burn in the presence of oxygen (air contains approximately 20% oxygen) to produce water and carbon dioxide. For example: methane (natural gas), which burns in air to produce water, carbon dioxide and maybe the most important additional product, heat. This energy in the form of heat is used to warm houses:



Another important example is the metabolism of glycogen. The combustion of saccharine (sugar) also creates carbon dioxide, water and energy:



This energy maintains the body temperature and is also used for the other energy requirements of the body.

It is common to think that combustion ends when all the oxygen is consumed. In this experiment, we will investigate it by burning a candle in a closed container while monitoring the oxygen content.

## Procedure

### Experiment setup



1. Remove the rubber protection from the Oxygen Sensor's cap. Unscrew the cap, fill half of it with the included liquid and then close it.


### Sensor setup

2. Connect the USB Module  to the USB port on your computer
3. Connect the oxygen sensor  to the USB Module.

#### Note:

**You must wait at least 5 minutes after connecting the sensor to the USB Module before beginning calibration and measurements.**

4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
5. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
6. This sets up the experiment parameters as follows:  
Experiment duration to 2 minutes  
Sampling rate to 10 per second  
Sensor mode: % in air.
7. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
8. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 2 minutes  
Set the sampling rate to 10 per second on the drop down menu.
9. Close the dialogue box.
10. Click the sensor's 'Setup'  icon in the sensors module box on the left side of the screen to open a dialogue box. Verify the mode is "% in air".
11. Whether you used the 'Load Configuration' function or not, click on the 'Calibration' icon . This will calibrate the sensor to a value of 20.9%.



12. Close the module setup dialogue box.

### Testing and measurements

13. Take the base with the side hole and insert the oxygen sensor in the hole, as in the picture. Put the candle in the center of the base and light it.



14. Click on the 'Run Experiment' icon again.

15. Do not stop the measurement; wait about 10 seconds and cover the candle (and the sensor) with the glass container.

16. Observe the changes in the oxygen percentage.

### Summary Questions

1. Explain why the candle is extinguished.
2. Observe the graph and determine the oxygen percentage when the candle is extinguished.
3. What combustible was used in this experiment?

### Challenge

1. Investigate which chemical products are obtained as a result of combustion reactions.
2. Make a table with the physical and chemical properties of oxygen.
3. Analyze and conclude.
4. You can repeat the experiment with glasses of different sizes to see how this affects the results.

# Endothermic and Exothermic Reactions 1

## Objective

- To study the enthalpy of reaction.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 

NOTE: Individual **NeuLog™** sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Materials

- Water
- Vinegar ( $\text{CH}_3\text{COOH}$ )
- Sodium bicarbonate  $\text{NaHCO}_3$
- A paper cup

## Discussion

In some reactions, the energy that must be absorbed to break the bonds in the reactants is less than the total energy released when new bonds are formed. This means that in the overall reaction, energy is released as either heat or light. This type of reaction is called an exothermic reaction. Another way of describing an exothermic reaction is that it is one in which the energy of the product is less than the energy of the reactants because energy has been released during the reaction.

In other reactions, the energy that must be absorbed to break the bonds in the reactants is more than the total energy that is released when new bonds are formed. This means that in the overall reaction, energy must be absorbed from the surroundings. This type of reaction is known as an endothermic reaction. Another way of describing an endothermic reaction is that it is one in which the energy of the product is greater than the energy of the reactants because energy has been absorbed during the reaction.

The difference in energy (E) between the reactants and the products is known as the **heat of the reaction**. It is also sometimes referred to as the **enthalpy** change of the system.

In this experiment, you will study the reaction of baking soda and vinegar:



## Procedure

### Experiment setup

1. Assemble a system as in the picture below.



2. Pour 50 ml vinegar into a cup and put it in front of the stand.

### Sensor setup

3. Connect the USB Module  to the USB port on your computer.
4. Connect the Temperature Sensor  to the USB Module using a Data Cable.
5. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
6. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
7. This sets up the experiment parameters as follows:  
 Experiment duration to 2 minutes  
 Sampling rate to 20 per second

8. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
9. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 2 minutes  
Set the sampling rate to 20 per second on the drop down menu.
10. Close the dialogue box.

### Testing and measurements

11. Insert the Temperature Sensor into the cup.
12. Click the 'Run Experiment' icon  to start the measurement.
13. Allow the temperature to stabilize.
14. Without stopping the measurement, slowly add 5 grams (approximately a teaspoon) of sodium bicarbonate.
15. Observe the changes in temperature.

### Summary Questions

1. What happens to the temperature during the reaction?
2. Is this an exothermic or an endothermic reaction?
3. Calculate the temperature change during the reaction.
4. Calculate the amount of heat energy lost in the reaction using the equation:  
 $q = m \times C_p \times \Delta T$ . The value of  $C_p$  is  $4.18 \text{ J/g}^\circ\text{C}$ .

### Challenge

1. Do you think it would be better to use a glass beaker instead of the paper cup in this experiment? Explain your answer.

## Endothermic and Exothermic Reactions 2

### Objective

- To study the enthalpy of reaction.
- To follow temperature changes during the dissolution of sodium hydroxide in water.

### Modules and Sensors

- Computer with NeuLog™ Software
- 2 Data Cables, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 
- pH Sensor: S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

### Materials

- Water
- 1 teaspoon solid Sodium hydroxide, about 2 grams
- A paper cup (a coffee cup with cover)

### Discussion

Many chemical reactions involve either the release or absorption of energy. In some reactions, the energy that must be absorbed to break the bonds in the reactants is less than the total energy released when new bonds are formed. This means that in the overall reaction, energy is released as either heat or light. This is called an exothermic reaction. Another way of describing an exothermic reaction is that it is one in which the energy of the product is less than the energy of the reactants because energy was released during the reaction.

In other cases, the energy that must be absorbed to break the bonds in the reactants is more than the total energy released when new bonds are formed. This means that in the overall reaction, energy must be absorbed from the surroundings. This is known as an endothermic reaction. Another way of describing an endothermic reaction is that it is one in which the energy of the product is greater than the energy of the reactants, because energy has been absorbed during the reaction.

The difference in energy (E) between the reactants and the products is known as **heat of the reaction**. It is also sometimes referred to as the **enthalpy** change of the system.


In this experiment, you will study the dissolution of sodium hydroxide in water and decide whether it is an exothermic or endothermic reaction.

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble a system as in the picture below.




2. Take the plastic cover of a coffee cup and make two holes one for the pH and the other for the Temperature Sensor (see picture).
3. Pour 100 ml water into the cup, put it in front of the stand and put the cover on the cup.
4. Prepare 1 teaspoon solid Sodium hydroxide, about 2 grams.

### Sensor setup

5. Connect the USB Module  to the USB port on your computer.
6. Connect the Temperature  and the pH  Sensors to the USB Module using the Data Cables.

7. Unscrew the cap of the storage solution of the pH Sensor and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
8. Introduce the pH Sensor  into the cup of water through the hole in the cover and stir carefully to remove the soaking solution.
9. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
10. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
11. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
12. This sets up the experiment parameters as follows:  
Experiment duration to 15 minutes  
Sampling rate to 60 per minute
13. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
14. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 15 minutes  
Set the sampling rate to 60 per minute on the drop down menu.
15. Close the dialogue box.
16. Whether you used the Load Configuration function or not, click on the 'Module Setup' icon  in the pH Sensor module box to open a dialogue box. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).



17. Insert the Temperature Sensor into the cup, through the hole in the cover.

### Testing and measurements

18. Click the 'Run Experiment' icon  to start the measurement.
19. Allow the temperature and the pH readings to stabilize.
20. Without stopping the measurement, open the cover of the cup and slowly add 1 teaspoon of solid Sodium hydroxide.
21. Observe the changes of temperature and pH.

### Summary Questions

1. What happens to the temperature during this reaction?
2. Is this an exothermic or endothermic reaction?
3. Calculate the temperature change during the reaction. You could use the 'Show Cursors'  icon. The values, at the bottom left-hand corner of the graph, indicate the X and Y for each cursor and the difference between X1 and X2 and between Y1 and Y2.
4. Calculate the pH change during the reaction (for a better resolution of the change in pH you can use the 'Show Cursors'  icon).
5. Did the pH change fast? Did the temperature change as fast as the pH?

### Challenge

1. Calculate the amount of heat energy released in the reaction using the equation:  
 $q = m \times C_p \times \Delta T$ . The value of  $C_p$  of water is  $4.18 \text{ J/g}^\circ\text{C}$ ,  $m$  is 100 g.
2. Try to dissolve different amounts of NaOH in water, calculate the heat of reaction in each case.

# Distillation – Part I

## Objective

To investigate the technique of distillation by separating the components of a mixture.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Distillation Kit: S02112
- Stand with utility clamp

## Materials

- Boiling stones
- Food coloring
- Ethanol

## Discussion

Distillation is a method for separating mixtures based on differences in their volatilities in a boiling liquid mixture. Distillation is a unit operation or a physical separation process, and not a chemical reaction.

Commercially, distillation has a number of uses such as the separation of crude oil into more fractions for specific uses such as transport, power generation and heating. Water is distilled to remove impurities, such as salt from seawater. Air is distilled to separate its components – notably oxygen, nitrogen, and argon – for industrial use. Distillation of fermented solutions has been used since ancient times to produce distilled beverages with higher alcohol content. The

premises where distillation is carried out, especially distillation of alcohol are known as a distillery.


The process of separating fractions of a mixture, boiling it at determined temperature intervals to obtain a quite defined composition of the product is called fractional distillation.

## Procedure

### Experiment setup

**Warning: You will be working with an open flame. Make sure not to have any flammable material in the proximity of the experiment. Be very careful when working with high temperatures and be sure to utilize proper laboratory safety techniques.**

1. Assemble the device as shown in the picture. This experiment does not require cooling with water.
  - a. In a distillation flask, pour 10 ml of Ethanol, some boiling stones and a pinch of coloring.
  - b. Put the distillation adaptor on the flask, screw the black cap tightly.
  - c. Put the condenser on the side tube of the adaptor and screw the black cap on the adaptor tightly.
  - d. Put the bended tube on the other side of the condenser and screw the black cap tightly.
  - e. Insert the Temperature Sensor in the distillation adaptor through the hole in the top cap up to where the distillation adaptor and the condenser meet (see picture). Screw this cap tightly to hold the sensor in place.
  - f. Put the assembled system in the stand, holding it through the utility clamp.
  - g. Unscrew the cap of the burner; fill the burner with ethanol and screw the cap back in place tightly.
  - h. The wick should be down so that it gets soaked with ethanol. After a few minutes, bring it up by turning the screw on the side of the cap. Feel the wick to make sure it is moist.



### Sensor setup

2. Connect the USB Module  to the USB port on your computer.
3. Connect the Temperature Sensor  to the USB Module with a Data Cable.
4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
5. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
6. This sets up the experiment parameters as follows:  
Experiment duration to 10 minutes  
Sampling rate to 60 per minute
7. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
8. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 10 minutes  
Set the sampling rate to 60 per minute on the drop down menu.
9. Close the dialogue box.

### Testing and measurements

10. Click on the 'Run Experiment' icon  to start the measurement.
11. Keep the flask about 5 cm above the flame. Slowly heat the mixture.
12. Turn off the burner at the end of the measurement (if all Ethanol evaporates before the time is over, stop the experiment).
13. At the end of the experiment, take apart the distillation system you assembled, wash all the parts with tap water and allow them to dry well before putting the distillation kit away.

## Summary Questions

1. Observe the graph and determine the evaporation point of Ethanol.
2. What was the color of the Ethanol that dripped into the beaker, was it colorless or did it have its original color?
3. Why?

## Challenge

1. Explain the different changes of state (evaporation, condensation, solidification, fusion and sublimation).
2. Give three examples of mixtures that could be separated through distillation.
3. Compare your different observations. Analyze and conclude.

## Distillation – Part 2

### Objective

- To investigate the technique of distillation by separating the components of a mixture.

### Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

### Equipment and Accessories

- Distillation Kit: S02112
- Stand with utility clamp

### Materials

- Boiling stones
- Food coloring
- Ethanol

### Discussion

Distillation is a method of separating mixtures based on differences in their volatilities in a boiling liquid mixture. Distillation is a unit operation or a physical separation process, and not a chemical reaction.

Commercially, distillation has a number of uses such as the separation of crude oil into more fractions for specific uses such as transport, power generation and heating. Water is distilled to remove impurities, such as salt from seawater. Air is distilled to separate its components – notably oxygen, nitrogen, and argon – for industrial use. Distillation of fermented solutions has been used since ancient times to produce distilled beverages with higher alcohol content. The premises where distillation is carried out, especially distillation of alcohol are known as a


distillery. The distillation that separates fractions of a mixture, boiling at determined temperature intervals to obtain a quite defined composition of the product is called fractional distillation.

## Procedure

### Experiment setup

**Warning: You will be working with an open flame. Make sure not to have any flammable material in the proximity of the experiment. Be very careful when working with high temperatures and be sure to utilize proper laboratory safety techniques.**

1. Assemble the device as shown in the picture. This experiment does not require cooling with water.
  - a. In a distillation flask pour about 5 ml of water with 5 ml of Ethanol, some boiling stones and a pinch of coloring.
  - b. Put the distillation adaptor on the flask, screw the black cap tightly.
  - c. Put the condenser on the side tube of the adaptor and screw the black cap on the adaptor tightly.
  - d. Put the bended tube on the other side of the condenser and screw the black cap tightly.
  - e. Insert the temperature sensor in the distillation adaptor through the hole in the top cap up to where the distillation adaptor and the condenser meet (see figure). Screw this cap tightly to hold the sensor in place.
  - f. Put the assembled system in the stand, holding it through the utility clamp.
  - g. Unscrew the cap of the burner; fill the burner with ethanol and screw the cap back in place tightly.
  - h. The wick should be down so that it gets soaked with ethanol. After a few minutes, bring it up by turning the screw on the side of the cap. Feel the wick to make sure it is moist.



### Sensor setup

2. Insert the USB Module  to the USB port on your computer.
3. Connect the Temperature Sensor  to the USB Module with a Data Cable.
4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
5. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
6. This sets up the experiment parameters as follows:  
 Experiment duration to 10 minutes  
 Sampling rate to 60 per minute
7. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
8. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
 Set the Experiment Duration to 10 minutes  
 Set the sampling rate to 60 per minute on the drop down menu.
9. Close the dialogue box.

### Testing and measurements

10. Click on the 'Run Experiment' icon  to start the measurement.
11. Keep the flask about 5 cm above the flame. Slowly heat the mixture.
12. Turn off the burner at the end of the measurement (if all the water evaporates before the time is over, stop the experiment).
13. At the end of the experiment, take apart the distillation system you assembled, wash all the parts with tap water and allow them to dry well before putting the distillation kit away.

## Summary Questions

1. Observe the graph and compare it with the one you received in Distillation – Part I. Describe the differences and explain them.
2. What compound was liberated first? Why?

## Challenge

1. Investigate the boiling point of Ethanol.
2. Try a distillation of a different mixture: 5 ml water with 5 ml Acetone.

# Boyle's Law

## Objective

- To measure the change in pressure as a function of change in temperature in a gas sample.

## Modules and Sensors

- Computer with NeuLog™ Software
- 2 Data Cables, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Temperature Sensor: S98242-31ND 
- Pressure Sensor: S98242-13ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Gas Law Device Kit: S02107  
(Volume of gas in the device is 0.055 L)
- Cables

## Discussion

An ideal gas is defined as one in which all collisions between atoms or molecules are perfectly elastic and in which there are no intermolecular attractive forces. One can visualize it as a collection of perfectly hard spheres which collide but which otherwise do not interact with each other. In such a gas any change in internal energy is accompanied by a change in temperature.

An ideal gas can be characterized by three parameters: volume, pressure and temperature for a certain amount of gas. The law that presents the relationship between volume, pressure and temperature is called the Ideal gas Law:  $PV = nRT$  where  $nR$  is constant when mass of the gas is constant; the temperature is in Kelvin degrees.

When temperature increases, the volume is constant, the gas expands and therefore the pressure increases. It is not necessary to know the exact value of the constant  $nR$  in order to


use the law; if the amount of gas and the volume are kept constant the relationship:  $P_1 / T_1 = P_2 / T_2$  must be fulfilled.

In this activity, we will change the temperature of the gas in a container while keeping the volume constant, in order to measure how it affects the pressure of the gas.

## Procedure

### Experiment setup

**Warning:** You will be working with an open flame. Make sure not to have flammable materials in the proximity of the experiment. Be very careful when working with high temperatures and be sure to utilize proper laboratory safety techniques.



Picture of experiment setup

1. Unpack Gas Law Device Kit (shown in the picture).
2. Put the candles under the device to heat up the gas in the chamber.

### Sensor setup

3. Connect the USB Module  to the USB port on your computer.
4. Connect the Pressure  and the Temperature  Sensors to the USB-100 module.
5. Insert the Pressure Sensor's inlet connector into one of the kit's perforated caps, as shown in the picture.
6. Insert the Temperature Sensor into one of the kit's perforated caps. Make sure that the third cap is sealed.
7. Make sure that the third cap is sealed and verify that there no air is leaking by pressing the syringe's plunger and releasing it so that it goes back up. If necessary, take out the syringe, bring the plunger to the 60 ml mark and put it back in place.
8. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
9. This sets up the experiment parameters as follows:
  - Sensor's units: Atm
  - A graph with the appropriate scale
  - A table
10. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
11. Click on the 'Experiment Setup' icon . This opens a dialogue box. Click the box next to table.
12. Close the dialogue box.
13. Click the 'Module Setup' icon  in the Pressure Sensor's module box to open a dialogue box.
  - Set the sensor's units by selecting the radio button next to Atm.
  - Set the graphs scale as follows:
  - Select the graph tabulator in the dialogue box.
  - Use the arrows to set the Y max. to 2 and the Y min. to 0.
  - The Y-axis position should be 1.

14. Click the 'Module Setup' icon  in the Temperature Sensor's module box to open another dialogue box.

Set the graphs scale as follows:

Select the graph tabulator in the dialogue box.

Use the arrows to set the Y max. to 70 and the Y min. to 10.

The Y-axis position should be -7.

### **Testing and measurements**

15. Record the temperature before starting the experiment. For the first measurement, click the 'Single Step' icon .

16. Start heating the gas using 3 candles. You can monitor the temperature increase by observing the temperature readings in the sensor's module box on the left side of the screen.



17. After every increase of 3 degrees in temperature, click the 'Single Step' icon  for a new measurement.

18. Continue the measurements until you reach 50 °C.

19. Blow the candles, remove them and stop the experiment. Do not to touch the container while it is still hot.

## Summary Questions

1. Draw a graph of pressure as a function of temperature. Click the Experiment Parameters icon  and select the tabulator Graph XY. Select the sensor radio button and then temperature from the drop down menu.
2. In order to calculate the linear fit for the points we received in the graph, click the 'Functions' icon .
3. Verify that the functions tabulator is active and also that the function which appears in the second drop down menu is linear fit. Click the 'Calculate Function' icon  so that you can see in the graph the corresponding line and in the functions window, the calculated linear function.

## Challenge

1. Based on the data and graph obtained for this experiment, express in words the relationship between gas pressure and temperature.
2. Explain this relationship using the concepts of molecular velocity and collisions of molecules.

# Producing Electricity

## Objective

- To determine the electrode combination that produces the best potential difference in a simple electrochemical cell.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- Voltage Sensor: S98242-7ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Red and black cables with crocodile connectors

## Materials

- Lemon or vinegar
- Electrodes made of available substances:
  - Carbon from a pencil (C)
  - Aluminum from aluminum foil (Al)
  - Iron from a nail (Fe)
  - Cooper from an electric cable (Cu)

## Discussion


In 1800 Alessandro Volta discovered that electricity could be produced by using small sheets of copper, zinc and cloth spacers soaked in an acid solution. He built a battery - the first apparatus capable of producing electricity.

Since then, batteries have been developed and have a great role in our daily life. A battery is a device that converts chemical energy directly to electrical energy and it is made of electrochemical cells. These often consist of two different materials inside a solution called electrolyte and are connected amongst them by a conductor.

In this experiment, you will build simple electrochemical cells using lemon or citric acid as the electrolyte and different materials as electrodes. You will measure the voltage produced when you introduce the different electrodes into the lemon and observe how some solids conduct electricity better than others.

## Procedure

### Experiment setup




1. Build a system as shown in the picture above.
2. Connect the crocodile connector of the red cable to the red cable of the Voltage Sensor (see picture).
3. Connect the crocodile connector of the black cable to the black cable of the Voltage Sensor (see picture).
4. Make sure the red and black connections do not touch each other (see picture).

The electrodes are: C, Fe, Cu and Al.

5. Make two separate cuts in the lemon (approximately 2 cm apart). If a lemon is not available, use vinegar in a cup instead.

### Sensor setup

6. Connect the USB Module  to the USB port on your computer.
7. Connect the Voltage Sensor  to the USB Module with a Data Cable.
8. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
9. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
10. This sets up a table.
11. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
12. Click the 'Experiment Setup' icon . This opens a dialogue box. Select the box next to Table.
13. Close the dialogue box.

### Testing and measurements

14. Insert the pencil (C) into one cut in the lemon and the nail (Fe) in the second cut. Connect the free red crocodile connector to the pencil and the free black crocodile connector to the nail.
15. When the reading stabilizes, click the 'Single Step' icon  to record the measurement in the table.
16. Write the electrode name for the voltage measurement in your table.
17. Remove the nail and insert the aluminum foil (Al). Connect the black crocodile to the aluminum foil. Click again the 'Single Step' icon  to record the new measurement in the table and in the graph.

18. Remove the aluminum foil and insert the electric cable (Cu). Connect the black crocodile to the electric cable. Click again the 'Single Step' icon  to record the new measurement in the table and in the graph.

Data Table:

Crocodile Conector		Average Voltage	Observations
Red end	Black end		
C	Fe		
C	Al		
C	Cu		

## Summary Questions

1. What combination of electrodes gives a higher voltage?
2. According to the previous answer and your own criteria, what combination would give the best electrochemical cell? Why?
3. Insert the best pair of electrodes again into the lemon, connect them to the crocodile connector and leave them for a while. Observe the voltage as time passes.
4. Write down the chemical equations that describe the voltage generation (oxidation-reduction processes) for each pair of electrodes.

## Challenge

1. Inverse the terminals and repeat the experiment. Conclude.

## Solar Cells

### Objective

- To study how a solar cell works.
- To investigate the correlation between the wavelength of the light and the photovoltaic voltage output.

### Modules and Sensors

- Computer with NeuLog™ Software
- 3 Data Cables, 10cm – S98242-48ND
- USB Module - S98242-45ND 
- Voltage Sensor - S98242-7ND 
- Light Sensor - S98242-20ND 
- Battery Module - S98242-44ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

### Equipment and Accessories

- A photovoltaic cell, (PV cell)
- Filters of different colors: red, yellow and dark blue

### Discussion

**Photovoltaic Cells** (or **PV**) are the field of technology and research related to the application of solar cells for energy by converting solar energy directly into electricity. Some materials exhibit a property known as the photoelectric effect which causes them to absorb photons of light and release electrons. These free electrons generate an electric current. Solar cells produce direct current (DC) electricity; an inverter can be used to change it to alternating current (AC) electricity. This electricity can be stored in batteries or other storage mechanisms to be used in absence of solar energy.

The photoelectric effect was first noted by a French physicist, Edmund Becquerel, in 1839. He found that certain materials produced small amounts of electric current when exposed to light. In 1905, Albert Einstein described the nature of light and the photoelectric effect on which photovoltaic technology is based, for which he later won a Nobel Prize in physics. Practical

photovoltaic (PV) cells were developed in 1954 when they were used to power toys. In 1958 they found wide acceptance as part of the space program after initial success on the Vanguard I satellite.

Solar cells are made from the same kind of semiconductor materials, such as silicon, used in the microelectronics industry. For solar cells, a thin semiconductor wafer is specially treated to form an electric field, positive on one side and negative on the other. When light energy strikes the solar cell, electrons are knocked loose from the atoms in the semiconductor material. If electrical conductors are attached to the positive and negative sides, forming an electrical circuit, the electrons can be captured in the form of an electric current -- that is, electricity. This electricity can then be used to power a load, such as a light or a tool.


In the first part of this experiment, we will measure the voltage produced by a photovoltaic cell when exposed to a table lamp and to solar radiation. The power of the solar cell is directly proportional to the voltage and the current.

In the second part, we will investigate the correlation between the wavelength of the light and the photovoltaic voltage output, using filters of different colors.

## Procedure


### Experiment setup

1. Assemble the experiment setup as shown in the picture.




2. Connect the red lead of the Voltage Sensor  to the red connector of the PV cell and the black lead of the sensor to the black connector of the PV cell.
3. To carry out the experiment indoors, use a table lamp to activate the PV cell. Place the Light Sensor  in front of the source of light.

**Sensor setup for indoors measurements:**

4. Connect the USB Module  to the USB port to your computer.
5. Using a Data Cables, connect the Voltage  Sensor to the USB Module and then the Light  Sensor to the Voltage Sensor, creating a chain.
6. Run the NeuLog™ software and check that the sensors are identified. If the software is already running, click the ‘Search for Sensors’ icon .
7. If this experiment was opened through the ‘Load Activity’ icon  located in the ‘Tools’ sub icon bar , click on the ‘Load Configuration’ icon .
8. This sets up the experiment parameters as follows:  
Experiment duration to 10 seconds  
Sampling rate to 10 per second  
Light sensor's range: 6000 Lux.
9. If you did not use the ‘Load Configuration’ function , click the ‘On Line Experiment’ icon  in the main icon bar.
10. Click on the ‘Experiment Setup’ icon . This opens a dialogue box:  
Set the Experiment Duration to 10 seconds  
Set the sampling rate to 10 per second on the drop down menu.
11. Close the dialogue box.

12. Click on the Light Sensor's 'Module Setup' icon  to open a dialogue window. Set the sensor's range by selecting the radio button next to 6,000 Lux.



13. Close the module setup dialogue box.


### Testing and measurements


14. Place the table lamp over the PV cell.
15. Hold the light sensor  in front of the source of light.
16. Click the run activity  icon to start the measurement for both sensors, light and voltage.
17. Click the 'Show Functions'  icon and select the Statistics tab. Select voltage with the **down-arrow**  of drop down menu and then click the 'Calculate Statistics'  icon to perform the calculation of the average.
18. Write the result in the indoors data table. Repeat for light average.
19. Repeat steps 14-16 for two more times and compare results.

### **Effect of the wavelength**

20. Put one colored filter over the light sensor, tilt the sensor towards the lamp and click the 'Run Experiment'  icon to start the measurement.
21. Average the light intensity as in step 17 and record the value in the data table.
22. Place the same filter over the PV cell and press the 'Run Experiment'  icon again to start a new measurement.
23. Average the voltage as in step 17 and record the value in the data table.
24. Apply steps 19-22 to each colored filter you have.

**Sensor setup for outdoors measurements:**

25. Connect the USB Module  to the USB port on your computer
26. Connect the Voltage  and the Light  Sensors to the USB module using Data Cables.
27. Since this activity uses Sun radiation, it must be performed outside. You will work in the Off Line mode of the NeuLog software to program the measurement conditions for the sensors to work disconnected from the computer (outside).
28. Click the Off-line experiment  icon in NeuLog's main icon bar.
29. Click the 'Module Setup' icon  in the sensors module box to open a dialogue box:  
Set the experiment duration to 10 seconds  
Set the light sensor's range by selecting the radio button next to 150,000 Lux



This procedure downloads the measurement conditions to the sensor and these are stored in its memory.

30. Close the module setup dialogue boxes.
31. Disconnect the sensors from the USB Module.
32. Connect the sensors to the Battery Module .

**Testing and measurements**

33. Tilt the PV cell towards the sun. Hold the Light Sensor at the same angle.
34. Collect the data by pressing the sensors' start/stop buttons.

**Note:**

You should see the sensors' red LED On during the measurement. When the LED turns Off, it means the experiment time is over. The measured data is stored in the sensors' memory.


35. Repeat steps 32- 33 for two more times and compare your results. You can save up to 5 experiments in the sensor's memory when using the Off Line Experiment mode.
36. At the end of the experiment outside, create a chain with the two sensors and connect them to the USB Module to down load the data into the computer as follows:

37. Click the search for sensors icon , then click the 'Off Line Experiment' icon in NeuLog's main icon bar. 

38. Click on the 'Load data from Sensors' icon  in the sub-icons bar to reveal a box with options:

**All (last experiments)**  
**Light 1**  
**Voltage 1**  
**Experiments**

39. Select **Experiments**. This will open the window: "Upload experiment number". See following picture:



40. Upload the measurements to your computer. The measurements were saved in the sensor's memory. Select the two radio buttons next to the number 5 (one for each sensor). This way, you are selecting the first measurements to be uploaded.

41. Click on the 'Load Data from Sensors' icon  inside the "upload experiment number" window to load the selected experiment. Observe the graphs.

42. For each measurement and each sensor, average the data as in step 17 and record the value in the outdoors data table.

**Effect of the wavelength**

43. Put a colored filter over the light sensor and tilt the sensor towards the Sun.
44. Collect data by pressing the Light Sensor's start/stop button.
45. Place the same filter over the PV cell and tilt the cell towards the Sun. Collect data by pressing the Voltage Sensors' start/stop button.
46. Apply steps 38-40 to each colored filter you have.
47. At the end of the experiments outside, connect the sensors back to the USB Module and download the data as in steps 35-36.
48. For each measurement and each sensor, average the data as in step 17 and record the value in the outdoors data table.

**Data Collection:**

**Indoors PV Voltage output**

	Voltage	Illumination (Lux)
1		
2		
3		
Average		

**PV Effect of the wavelength**

	Voltage	Illumination (Lux)
No Filter		
Yellow filter		
Red filter		
Blue Filter		

**Outdoors PV Voltage output**

	Voltage	Illumination (Lux)
1		
2		
3		
Average		

**PV Effect of the wavelength**

	Voltage	Illumination (Lux)
No Filter		
Yellow filter		
Red filter		
Blue Filter		

**Summary Questions**

1. Calculate the average voltage and illumination values for indoors and outdoors experiments.
2. What is the difference in voltage output for each of the trials for each of the filters?
3. Which of the colors had the highest power output? Which had the lowest output?

**Challenge**

1. Determine the effect of the amount of light shining on a PV cell on the voltage output.
2. What can you conclude about the effect of wavelength on power output?

# Buffered Vitamin C Tablets

## Objective

- To study a buffered solution by monitoring the pH value of two different types of vitamin C tablets.
- Compare a buffered vitamin C with a regular vitamin C tablet.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm – S98242-48ND
- USB Module - S98242-45ND 
- pH Sensor - S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Ring stand with utility clamp

## Materials

- Two different vitamin C 500 mg tablets, one regular and one buffered
- Two paper cups
- Distilled or tap water


## Discussion

In chemistry, a buffer is a system, usually an aqueous (water) solution, which resists having its pH changed when an acid or a base is added to it.

Buffer solutions are used as a means of keeping the pH at a nearly constant value in a wide variety of chemical applications. Many life forms thrive only in a relatively small pH range. An example of a buffer solution is blood.

There are two common kinds of buffer solutions: those that contain a weak acid plus one of its salts (e.g., acetic acid plus sodium acetate) and solutions that contain a weak base plus one of its salts (e.g., ammonia plus ammonium chloride).

There is equilibrium in a solution between a weak acid, HA, and its conjugate base, A<sup>-</sup> :



- When hydrogen ions (H<sup>+</sup>) are added to the solution, the equilibrium moves to the left, as there are hydrogen ions (H<sup>+</sup> or H<sub>3</sub>O<sup>+</sup>) on the right-hand side of the equilibrium expression.
- When hydroxide ions (OH<sup>-</sup>) are added to the solution, the equilibrium moves to the right, as hydrogen ions are removed in the reaction (H<sup>+</sup> + OH<sup>-</sup> → H<sub>2</sub>O).

Thus, in both cases, some of the added reagent is consumed in shifting the equilibrium in accordance with Le Chatelier's principle and the pH changes by less than it would if the solution was not buffered.


In this experiment we will use a pH Sensor to follow the behavior of regular vitamin C tablets as compared to that of a buffered vitamin.

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.


1. Assemble a system as the one shown in the picture below:



2. Pour 100 ml of distilled or tap water into a clean paper cup.
3. Put the cup in front of the stand.

### Sensor setup

4. Connect the USB Module  to the USB port on your computer.
5. Connect the pH Sensor  to the USB Module using a Data Cable.
6. Unscrew the cap of the storage solution and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
7. Put the pH Sensor in the tube holder on the stand.
8. Introduce the sensor into the water and stir carefully to remove the soaking solution.
9. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
10. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
11. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
12. This sets up the experiment parameters as follows:  
Experiment duration to 15 minutes  
Sampling rate to 30 per minute
13. If you did not use the 'Load Configuration' function , click the 'On Line Experiment' icon  in the main icon bar.
14. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 15 minutes  
Set the sampling rate to 30 per minute on the drop down menu.
15. Close the dialogue box.
16. Whether you used the 'Load Configuration' function or not, click on the 'Module Setup' icon  in the sensors module box to open a dialogue box.
17. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).



18. Close the module setup dialogue box.

**Testing and measurements**

- 19. Click the 'Run Experiment' icon  to start the measurement.
- 20. Allow the pH to stabilize. Record the initial pH value of distilled or tap water.
- 21. Without stopping the measurement, put the regular 500 mg vitamin C tablet. Record the pH value every 5 minutes.
- 22. At the end of the experiment, save the data by clicking the 'Freeze Current Graphs' icon .
- 23. Clean the pH Sensor with tap water.
- 24. Put 100 ml of tap water into a clean paper cup; add the buffered 500 mg vitamin C tablet.
- 25. Repeat steps 10-12, record the data.

**Data Collection:**

Time(min)	pH regular 500 mg vitamin C	pH Buffered 500 mg vitamin C

## Summary Questions

1. Print a graph for the pH vs. time for both vitamin C experiments.
2. Using the pH values in the data table, compare the change in pH ( $\Delta\text{pH}$ ) during time intervals. Which vitamin-C tablet had the largest decrease in pH each time? Explain.

## Challenge

1. Do you think buffered vitamins are better? Explain.

# Study of Slow-Release Food Supplements

## Objective

- To monitor how pH values change for two different types of vitamin C tablets dissolved in water over a ten minutes test.
- Study the benefit of slow-release food supplements using slow-release and regular vitamin C tablets dissolved in water.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm – S98242-48ND
- USB Module - S98242-45ND 
- pH Sensor - S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Materials

- Regular vitamin C tablet
- Slow release vitamin C tablet
- A cup
- Water
- **Note:** if possible, use effervescent vitamin C tablets.

## Discussion

When you take a regular vitamin pill, it is rapidly dissolved in the stomach and absorbed in the bloodstream. However, this rapid dissolution of the supplement increases the chances that your body will excrete significant traces of the beneficial vitamins. On the other hand, time-release vitamins use micro-technologies to dissolve much more slowly in the stomach. The common method for time-release is coating the vitamins, minerals or medicine with a polymer so they are released incrementally into the bloodstream and are absorbed over an extended period of time rather than just a few minutes. This increases the chance that a significant percentage of the vitamins or medicine will be absorbed in the body.

In this experiment we will use a pH sensor to follow an example of a slow-release food supplement by dissolving regular and slow release vitamin C tablets in water.

## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble a system as the one in the picture below.



2. Pour 100 ml of water into a cup.

### Sensor setup

3. Connect the USB Module  to the USB port on your computer.
4. Connect the pH Sensor  to the USB Module with a Data Cable.
5. Unscrew the cap of the storage solution and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
6. Introduce the sensor into the water and stir carefully to remove the soaking solution.
7. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .

8. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
9. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
10. This sets up the experiment parameters as follows:  
Experiment duration to 10 minutes  
Sampling rate to 30 per minute
11. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
12. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 10 minutes  
Set the sampling rate to 30 per minute on the drop down menu.
13. Close the dialogue box.
14. Whether you used the 'Load Configuration' function or not, click the 'Setup' icon  on the sensors module box to open a dialogue box. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used). Set the experiment duration to 10 minutes and the sampling rate to 30 per minute. Your screen should look as the following:



15. Close the module setup dialogue box.

### Testing and measurements

16. Click on the 'Run Experiment' icon  to start the measurement.
17. Without stopping the measurement, drop the regular vitamin C tablet into the cup. Swirl the cup from time to time.
18. Observe the graph of the measurement.
19. Click the 'Freeze Current Graph' icon .
20. Repeat steps 16 to 18 using the slow release vitamin C tablet.

### Summary Questions

1. Compare your graphs.
2. Compare the change in pH ( $\Delta\text{pH}$ ). What could you say about the rate of pH changes? For which kind of vitamin tablet the pH decreases faster?

### Challenge

1. Explain the benefit of slow release supplements to our body.

# Seed Respiration

## Objective

- To compare the production of carbon dioxide by germinating seeds and dormant seeds.
- To determine whether germinating seeds and non-germinating seeds respire and what factors can affect them.
- To determine the effect of temperature on seed germination.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- CO2 Sensor: S98242-3ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- A lamp
- A glass
- A plate
- Special glass bottle for the Carbon Dioxide sensor

## Materials

- 30 germinating seeds (lentils or any other seeds that are easy to germinate)
- 30 non-germinating seeds (lentils).

## Discussion

In some aspects, plants are not very different from people. Although plants can produce sugars from the Sun's energy, when plants need energy they have to metabolize their stored sugars through cellular respiration, just like we do.

A plant's respiration rate is not constant and depends on many factors. After a seed is separated from the plant, it usually goes into a resting period called dormancy. During

dormancy, the seed waits until conditions are just right before it begins the germination process.

During this time, seeds cannot produce their own food because they have no leaves. Therefore, in order for a seed to stay alive or grow it needs to use stored energy reserves and undergo cellular respiration. Have you ever wondered why seeds and nuts have so many calories? The seed will use these calories to survive during dormancy and to begin germination.

Cellular respiration occurs in every cell in both plants and animals and is essential for daily living. Cellular respiration is an exergonic reaction, which means it produces energy. It is also a catabolic process - it breaks down polymers into smaller, more manageable pieces. The ultimate goal of cellular respiration is to take carbohydrates and disassemble them into glucose molecules, and then use this glucose to produce energy-rich ATP molecules.

The general equation for cellular respiration is: (the energy released is about 36-38 molecules of ATP).



Seeds undergo cell respiration during germination. In this activity you will use the CO<sub>2</sub> gas sensor to monitor the carbon dioxide produced by seeds during cell respiration.

You will test both germinating and dormant seeds. In addition, you will test the effect of two different temperatures on cell respiration; other factors such as light and moisture will also be investigated.

## Procedure

### Experiment setup



1. Put about 30 non-germinating seeds, (lentils) into the special glass bottle.

### Sensor setup

2. Connect the USB Module  to the USB port on your computer.
  3. Connect the Carbon Dioxide sensor  to the USB Module using a Data Cable.
  4. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
  5. If this experiment was opened through the 'Load Activity' icon  located in the Tools sub icon bar , click the 'Load Configuration' icon .
  6. This sets up the experiment parameters as follows:  
Experiment duration to 10 minute  
Sampling rate to 50 per minute
  7. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
  8. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 10 minute  
Set the sampling rate to 50 per minute on the drop down menu.
  9. Close the dialogue box.
  10. Whether you used the Load Configuration function or not, click the sensor's 'Setup'  icon in the sensors module box on the left side of the screen to open a dialogue box and then click the 'Calibration' icon . This will calibrate the sensor to a value of 380 ppm. An alternative way of calibration is to press the sensor's calibration button continuously (3 seconds).
- Note:** Since our measurements are relative, the sensor can be calibrated in the classroom. For quantitative measurements, the calibration should be performed outside where there is fresh air.
11. Note that the measured amount of CO<sub>2</sub> appears in the sensor's module box. Close the module setup dialogue box.

### Testing and measurements

12. Introduce the sensor into the special glass bottle with the seeds.
13. Wait about 3 minutes for the reading to stabilize.

14. Click on the 'Run Experiment' icon  to start the measurement.
15. After the measurement, remove the sensor from the glass bottle and let it re-adjust for 3 minutes; observe how the CO<sub>2</sub> reading in the sensor's module box changes back to the initial value.
16. Save the data.  
**Note:** For a better resolution, click the 'Zoom Fit'  icon.
17. To determine the rate of respiration, click the 'Show Cursors' icon , move the two cursors, one to the part of the experiment where the increase starts to be linear and the second to the end of the measurement.
18. Click the 'Show Functions' icon , click the functions tab and select "Linear Fit (between cursors)" from the second drop down menu. Click the 'Calculate Function' icon  to perform the calculation; in the box you will see the formula with the best fit line. The slope of the graph is the rate of respiration. Record the result in your table.
19. Remove the non-germinating seeds from the bottle.

**Test germinating seeds:**

20. Obtain 30 germinating seeds and blot them dry between two pieces of paper towel; put them in the special glass bottle.
21. Introduce the sensor back in the glass bottle.
22. Repeat steps 13-19 for the germinating seeds.

**Test the effect of light:**

23. Introduce the sensor back into the special glass bottle. Fill a glass with water and place it between the lamp and the special bottle with the seeds and the sensor. (The glass with water protects the seeds from the heat of the lamp.)
24. Place the lamp as close to the glass with water as possible but without touching it. Note the time and turn on the lamp. The lamp should be on for 5 minutes before you begin the data collection.
25. Repeat steps 13-19 for the light irradiated germinating seeds.

**Test the effect of cold temperature:**

26. Introduce the sensor back into the special glass bottle.

27. Fill a deep plate with ice cold water and put the special bottle with the germinating seeds in it so that most of the bottle is in the water.
28. Repeat Steps 13–19 to collect data with the cold germinating seeds.

**Data Collection:**

<i>Seeds</i>	<i>CO<sub>2</sub> Rate of respiration (ppm/sec)</i>
<i>Dormant, room temperature</i>	
<i>Germinating, room temperature</i>	
<i>Germinating, room temperature with light</i>	
<i>Germinating, cool temperature</i>	

### Summary Questions

1. Did cell respiration occur in the seeds? Explain.
2. Compare the rate of CO<sub>2</sub> production of germinating and dormant seeds; what is the effect of germination on the rate of cell respiration?
3. How does temperature affect the CO<sub>2</sub> production rate for germinating seeds (room temperature compared to cold germinating)?
4. How does light affect the CO<sub>2</sub> production rate for germinating seeds at room temperature?

### Challenge

1. What other conditions beside temperature could affect the rate of respiration?
2. Grow your own seeds (lentils, beans, etc.). Measure the respiration rate at different stages of the germinating process.

# Acid Rain

## Objective

- To study the acid rain phenomenon.

## Modules and Sensors

- Computer with NeuLog™ Software
- Data Cable, 10cm: S98242-48ND
- USB Module: S98242-45ND 
- pH Sensor: S98242-24ND 

NOTE: Individual NeuLog™ sensors and connection hardware are sold individually or as part of economical Curriculum Starter Packs that include sensors and accessories in a convenient custom case.

## Equipment and Accessories

- Acid Rain Kit: S02157

## Materials

- Vinegar ( $\text{CH}_3\text{COOH}$ )
- Sodium bicarbonate ( $\text{NaHCO}_2$ )
- Water

## Discussion

Acid rain is a phenomenon associated with the development of urban and industrial areas. It consists of the incorporation of chemical compounds such as carbon dioxide, sulfur dioxide and nitrogen oxides to rain water. These compounds are emitted by cars' exhausts, factories and thermoelectric centrals.


When these oxides are in contact with humidity in the atmosphere, they are transformed into secondary polluting agents that form solutions of Carbon, Sulfuric and Nitric acid. Rain carries these compounds to the Earth's surface, deposits them on the soil and in water bodies.

The persistent fall of acid rain damages lakes, rivers and underground waters, causing the death of fish and other organisms in aqueous ecosystems. It acidifies and demineralizes soils, damages forests, national parks and reserves and causes low producing crops. Also, the

acidity of the water deteriorates archeological zones, historical monuments, buildings and metallic structures. The dimension of the damage depends on the degree of acidity.

In this experiment, you will produce one of the gases responsible of acid rain - CO<sub>2</sub> - and will quantify the variation in the pH induced in the water when it dissolves.

CO<sub>2</sub> and H<sub>2</sub>O (water) create acid H<sub>2</sub>CO<sub>3</sub> according to the following formula:




## Procedure

### Experiment setup

**Warning:** Please note that the bottom part of the pH Sensor consists of a fragile crystal sphere. Even though it has a plastic protection, be careful not to break it.

1. Assemble a system like the one shown below.



2. Put a cup with 100 ml water on the table.
3. Open the gas generator, add two spoonfuls of sodium bicarbonate and close the cap tightly.
4. Insert the hose attached to gas generator into the water.

### Sensor setup

5. Connect the USB Module  to the USB port on your computer.
6. Connect the pH Sensor  to the USB Module using a Data Cable.
7. Unscrew the cap of the storage solution and take out the sensor. Raise the cap to the top of the sensor. Put the solution aside.
8. Introduce the sensor into the water and stir carefully to remove the soaking solution.
9. Run the NeuLog™ software and check that the sensor is identified. If the software is already running, click the 'Search for Sensors' icon .
10. Observe the sensor's module box on the left side of the screen. Allow the pH value to stabilize.
11. If this experiment was opened through the 'Load Activity' icon  located in the 'Tools' sub icon bar , click on the 'Load Configuration' icon .
12. This sets up the experiment parameters as follows:  
Experiment duration to 5 minutes  
Sampling rate to 10 per second
13. If you did not use the Load Configuration function , click the 'On Line Experiment' icon  in the main icon bar.
14. Click on the 'Experiment Setup' icon . This opens a dialogue box:  
Set the Experiment Duration to 5 minutes  
Set the sampling rate to 10 per second on the drop down menu.
15. Close the dialogue box.
16. Whether you used the Load Configuration function or not, click on the 'Module Setup' icon  in the sensor module box to open a dialogue box.

17. Click on the 'Calibration' icon  to set the value at 7. (For a more accurate calibration, a buffer solution pH = 7 should be used).



18. Close the module setup dialogue box.

### Testing and measurements

19. Click on the 'Run Experiment' icon  to start the measurement.
20. Without stopping the measurement, slowly add about 50 ml of vinegar through the funnel. Swirl the gas generator from time to time.
21. You should see gas bubbles in the water.
22. Observe the pH changes in the graph on the computer screen.

## Summary Questions


1. We introduced gas into the water. What happens to the water, does it stay the same, becomes acid or basic according to the pH readings?
2. Does the pH level change during the experiment?

## Challenge

1. What are the main primary polluting agents that originate acid rain?
2. Develop the equation from the chemical reaction of acetic acid and sodium bicarbonate.



3. Develop the chemical reaction of acid rain from Sulfur trioxide plus water.



4. Develop the chemical reaction of the results of this activity.
5. What was the pH level obtained in the results?
6. Investigate the alterations to the environment produced by acid rain.
7. Analyze and conclude.