

Balance Buying Guide




MS-TS, ML-T, and ME-TE





METTLER TOLEDO recently released new basic and entry level balances: the MS-TS, ML-T, and ME-TE. Each boast a unique set of features and value, and together provide more solutions than ever before. Learn more about the families below, along with the previous generation models as reference to assist you in your selection.

Need Five Star Performance? ★ ★ ★ ★ ★

METTLER TOLEDO Excellence Balances (not shown here) can support the highest needs for safety, efficiency, quality assurance, and ease of compliance. Please visit fishersci.com/mettlertoledo or contact your local representative for details on anything you see here as well as our full lineup of Excellence balances and more advanced weighing solutions.

The New Generation															
MS-TS Balances					ML-T Balances					ME-TE Balances					
Sophistication Made Simple The MS-TS packs high-end features into an easy-to-use instrument. An extra-large 7" color touchscreen display, full-metal housing, and quick-lock draftshield make this a sophisticated balance that's dependable for the everyday user.					Big Features Go Compact The ML-T balance shares the solid performance of the MS-TS, but with a smaller touchscreen and footprint. Internal calibration, on-screen guidance and a metal bottom housing make it a strong and reliable addition to any lab.					For Your Essential Needs The ME-TE balance is an affordable option for those seeking the essentials. A touchscreen and built-in applications combined with a standard draft shield and external calibration make it a quality balance at a friendly price point.					
															
Performance	★	★	★	★	★	★	★	★	★	★	★	★	★	★	★
	★	★	★	★	★	★	★	★	★	★	★	★	★	★	★
	★	★	★	★	★	★	★	★	★	★	★	★	★	★	★
	★	★	★	★	★	★	★	★	★	★	★	★	★	★	★

The Previous Generation (Phased Out 2015)															
MS Balances								ML Balances							
Long-Lasting Quality Looking to replace your METTLER TOLEDO MS Balance? The new MS-TS and ML-T are the closest replacements to this popular product family.								Basic Weighing Looking to replace your METTLER TOLEDO ML Balance? For ML Analytical balances, the ML-T is your closest replacement. Precision (toploading) owners should check out either the ML-T or the ME-TE (external calibration only).							
															

For pricing and more info, visit fishersci.com/mettlertoledo

METTLER TOLEDO